

MET/CAL 10.1.0 HOT-FIX README

This document includes a list of updates included in the MET/CAL 10.1.0 Hot-Fix.

Note: This Hot-Fix only applies to MET/CAL version 10.1.0. Do not attempt to apply this update to any other version of MET/CAL!

This Hot-Fix applies to both MET/TEAM and MET/CONNECT systems. Refer to the table below for a list of issues that are resolved by this Hot-Fix.

APPLYING HOT-FIX TO MET/CAL VERSION 10.1.0

To apply this Hot-Fix, it will be necessary to shut down MET/CAL. This Hot-Fix should be applied to all MET/CAL workstations. This Hot-Fix only applies to MET/CAL Runtime application. The MET/CAL Editor application does not need to be updated.

On each MET/CAL workstation:

1. Make sure MET/CAL is not running.
2. Extract the files from the Hot-Fix archive.
3. Using File Explorer, navigate to the folder for the MET/CAL Runtime application (typically C:\Program Files (x86)\Fluke\METCAL).
4. Copy the **Metcal.Common.dll**, **MetConnect.Storage.dll**, and **Metteam.Storage.dll** files that were extracted from the Hot-Fix archive to the \METCAL folder, replacing the existing files.

Updating these files does not change the MET/CAL version. It will still appear as version 10.1.0.51 on the **About** dialog.

MET/CAL 10.1.0 HOT-FIX CHANGES

Key	Issue Type	Description
MC-1312	FIXED	Resolved an issue that caused an exception when writing results to the database when a procedure required a new Unit record to be created in the database.